

FGFR4-IN-23

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	FGFR4-IN-23 (compound 48c) is an FGFR4 inhibitor with an IC50 value of 2.9 nM for the FGFR4WT kinase. It effectively inhibits the phosphorylation of FGFR4, AKT, and MAPK.
Targets(IC50)	FGFR,Akt,AMPK

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481