

PPAR $\gamma$  modulator-3

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	PPAR $\gamma$ modulator-3 (Compound 11) is a modulator of the peroxisome proliferator-activated receptor $\gamma$ (PPAR $\gamma$ ) with a KD value of 186 nM. It holds potential for research into insulin resistance (IR)-related conditions, including type 2 diabetes mellitus (T2DM) and metabolic syndrome.
Targets(IC50)	PPAR

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481