

DSPE-PEG2000-TH

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	DSPE-PEG2000-TH is a PEG conjugate made from DSPE and a pH-responsive cell-penetrating peptide (TH). In acidic environments—such as the tumor microenvironment—TH becomes activated, allowing for the selective delivery of small molecules, oligonucleotides, and proteins into tumor cells.
Targets(IC50)	Liposome

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481