

DSPE-PEG5000-CCK8

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	DSPE-PEG5000-CCK8 is a PEGylated compound comprised of DSPE and Cholecystokinin-8 (CCK8). Cholecystokinin-8 exhibits peptide activity that regulates gallbladder contraction and the functioning of the digestive system. DSPE-PEG5000-CCK8 can be applied in drug delivery systems.
Targets(IC50)	Liposome

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481