

BTR2000

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Keep away from direct sunlight

Storage:

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	BTR2000 is a specific ligand for KLHL20, capable of binding with high affinity to the Kelch domain of KLHL20. It serves as a ligand for E3 ligase and is utilized in the synthesis of BTR2004.
Targets(IC50)	Ligands for E3 Ligase

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481