

DSPE-PEG1000-RVG29

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	DSPE-PEG1000-RVG29 is a PEG compound composed of DSPE and rabies virus glycoprotein 29 (RVG29). RVG29 can specifically bind to nicotinic acetylcholine receptors (nAChR) at the blood-brain barrier (BBB) and cross the BBB.
Targets(IC50)	AChR,Liposome

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481