

E3 Ligase Ligand-linker Conjugate 182

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: **Keep away from direct sunlight**
Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	E3LigaseLigand-linker Conjugate 182 is a type of E3LigaseLigand-Linker Conjugate composed of a CRBN ligand and a PROTAC linker. E3LigaseLigand-linker Conjugate 182 can be utilized in the design of PROTACs, such as [LZ-07].
Targets(IC50)	E3 Ligase Ligand-Linker Conjugates

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481