

BTR2000-Gly-Gly-Gly

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Keep away from direct sunlight  
Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
Actual storage temperature shall be subject to the COA.

## Biological Description

Description	BTR2000-Gly-Gly-Gly is a conjugate of an E3 ligase ligand (KLHL20) and a linker, which is utilized in the synthesis of PROTAC BTR2004. BTR2004 functions as a PROTAC degrader targeting BRD2/3/4.
Targets(IC50)	E3 Ligase Ligand-Linker Conjugates

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481