

LC3B ligand 2

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: **Keep away from direct sunlight**
Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	LC3B ligand 2 is a ligand for E3 ubiquitin ligase. It can be connected to a target protein ligand through a linker, forming PROTACPCSK9 autophagic degrader 2. PROTACs are capable of inducing the ubiquitination and degradation of oncogenic proteins.
Targets(IC50)	Ligands for E3 Ligase

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481