

DSPE-PEG1000-F3

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	DSPE-PEG1000-F3 is a PEG compound made from DSPE and a targeting peptide called nucleolin (F3). The F3 peptide can specifically bind to cell surface nucleolin, enabling efficient transport from the cell surface to the nucleus. DSPE-PEG1000-F3 is applicable in drug delivery.
Targets(IC50)	Liposome

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481