

DSPE-PEG2000-F3

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	DSPE-PEG2000-F3 is a PEG compound composed of DSPE and the nucleolin-targeting peptide (F3). The F3 peptide can specifically bind to nucleolin on the cell surface, facilitating efficient transport from the cell surface to the nucleus. DSPE-PEG2000-F3 is applicable for drug delivery.
Targets(IC50)	Liposome

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481