

## FGFR1 inhibitor-18

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	FGFR1 inhibitor-18 (compound 16) is a potent inhibitor of FGFR1 with an IC50 of 1.31 nM. It effectively occupies the ATP binding site of the target FGFR1.
Targets(IC50)	FGFR

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481