

## Goat Anti-Rat IgG H&L Secondary Antibody-AF594

### Product Details

Ig Type: IgG  
Reactivity: Rat  
Purification: Protein G purified

### Applications

Application: FCM,ICC/IF,IF  
Recommended IF: 1:200-1000; FCM: 1:50-200; ICC/IF: 1:100-1000

### Properties

Stability & Storage: Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.  
Shipping: Shipping with blue ice.

### Antigen Details

Immunogen: Native Rat IgG  
Synonyms: Goat Anti-Rat IgG H&L (AF594);Immunoglobulin G

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481