

## Donkey Anti-Rabbit IgG H&L Secondary Antibody-AF555



### Product Details




Ig Type:	IgG
Reactivity:	Rabbit
Purification:	Protein G purified

### Applications

#### Verified Activity:

1. Paraformaldehyde-fixed, paraffin embedded Rat Kidney; Antigen retrieval by boiling in sodium citrate buffer (pH6.0) for 15 min; The section was incubated with Vimentin Polyclonal Antibody, Unconjugated at 1:200 overnight at 4°C, followed by a conjugated Donkey Anti-Rabbit IgG antibody (TMAB-02007AF555) for 90 minutes, and DAPI for nuclei staining.
2. Paraformaldehyde-fixed, paraffin embedded Human Kidney; Antigen retrieval by boiling in sodium citrate buffer (pH6.0) for 15 min; The section was incubated with Vimentin Polyclonal Antibody, Unconjugated at 1:200 overnight at 4°C, followed by a conjugated Donkey Anti-Rabbit IgG antibody (TMAB-02007AF555) for 90 minutes, and DAPI for nuclei staining.
3. Paraformaldehyde-fixed, paraffin embedded Mouse Kidney; Antigen retrieval by boiling in sodium citrate buffer (pH6.0) for 15 min; The section was incubated with Vimentin Polyclonal Antibody, Unconjugated at 1:200 overnight at 4°C, followed by a conjugated Donkey Anti-Rabbit IgG antibody (TMAB-02007AF555) for 90 minutes, and DAPI for nuclei staining.
4. Paraformaldehyde-fixed, paraffin embedded Rat Small Intestine; Antigen retrieval by boiling in sodium citrate buffer (pH6.0) for 15 min; The section was incubated with Vimentin Polyclonal Antibody, Unconjugated at 1:200 overnight at 4°C, followed by a conjugated Donkey Anti-Rabbit IgG antibody (TMAB-02007AF555) for 90 minutes, and DAPI for nuclei staining.
5. Paraformaldehyde-fixed, paraffin embedded Mouse Small Intestine; Antigen retrieval by boiling in sodium citrate buffer (pH6.0) for 15 min; The section was incubated with Vimentin Polyclonal Antibody, Unconjugated at 1:200 overnight at 4°C, followed by a conjugated Donkey Anti-Rabbit IgG antibody (TMAB-02007AF555) for 90 minutes, and DAPI for nuclei staining.





Application: FCM,ICC/IF,IF

Recommended IF: 1:100-1000; ICC/IF: 1:100-1000; FCM: 1:100-1000

### Properties

Stability & Storage: Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.

Shipping: Shipping with blue ice.

### Antigen Details

Immunogen: Native Rabbit IgG

Synonyms: Immunoglobulin G;Donkey Anti-Rabbit IgG H&L (AF555)

### Research Background

Immunoglobulin G (IgG), is one of the most abundant proteins in serum with normal levels between 8-17 mg/mL in adult blood. IgG is important for our defence against microorganisms and the molecules are produced by B lymphocytes as a part of our adaptive immune response. The IgG molecule has two separate functions; to bind to the pathogen that elicited the response and to recruit other cells and molecules to destroy the antigen. The variability of the IgG pool is generated by somatic recombination and the number of specificities in an individual at a given time point is estimated to be 1011 variants.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481