

Goat Anti-Rabbit IgG H&L whole serum Secondary Antibody

Product Details

Ig Type:	IgG
Reactivity:	Rabbit
Purification:	Unpurified

Properties

Stability & Storage:	Store at -20°C or -80°C for 2 years (lyophilized powder). Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Native Rabbit IgG
Synonyms:	Rabbit Immunoglobulin G;Rb IgG

Research Background

immunoglobulin G; G1m marker; G2m marker; G3m marker; G4m marker; HDC; Heavy chain disease protein; Ig gamma 1 chain C region; Ig gamma 2 chain C region; Ig gamma 3 chain C region; Ig gamma 4 chain C region; IGHG 1; IGHG 2; IGHG 3; IGHG 4; IGHG1; IGHG2; IGHG3; IGHG4; Immunoglobulin heavy constant gamma 1; Immunoglobulin heavy constant gamma 2; Immunoglobulin heavy constant gamma 3; Immunoglobulin heavy constant gamma 4; MGC117419; MGC45809.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481