

## Goat Anti-Rabbit IgG H&L Secondary Antibody-SA

### Product Details

Ig Type:	IgG
Reactivity:	Rabbit
Purification:	Protein G purified

### Properties

Stability & Storage: Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.

Shipping: Shipping with blue ice.

### Antigen Details

Immunogen:	Native Rabbit IgG
Synonyms:	Rb IgG;Rabbit Immunoglobulin G

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481