

## Mouse Anti-Rabbit IgG H&L Secondary Antibody-AF594

### Product Details

Ig Type:	IgG
Reactivity:	Rabbit
Purification:	Protein G purified

### Applications

Application:	FCM,ICC/IF,IF
Recommended	IF: 1:100-1000; ICC/IF: 1:100-1000; FCM: 1:100-1000

### Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	Native Rabbit IgG
Synonyms:	Immunoglobulin G; Mouse Anti-Rabbit IgG H&L (AF594)

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481