

## Goat Anti-Mouse IgG H&L Secondary Antibody-SA

### Product Details

Ig Type:	IgG
Reactivity:	Mouse
Purification:	Protein G purified

### Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	Native Mouse IgG
Synonyms:	Mo IgG; Mouse Immunoglobulin G

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481