

Goat Anti-Human IgG H&L Secondary Antibody-AF594

Product Details

Ig Type: IgG
Reactivity: Human
Purification: Affinity purified by Protein G, nonspecific adsorbed

Applications

Application: FCM,ICC/IF,IF
Recommended IF: 1:100-500; FCM: 1:50-200; ICC/IF: 1:100-500

Properties

Stability & Storage: Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping: Shipping with blue ice.

Antigen Details

Immunogen: Native Human IgG
Synonyms: Goat Anti-Human IgG H&L (AF594);Immunoglobulin G

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481