

Rabbit Anti-Chicken IgG H&L whole serum Secondary Antibody

Product Details

Ig Type:	IgG
Reactivity:	Chicken
Purification:	Unpurified

Properties

Stability & Storage:	Store at -20°C or -80°C for 2 years (lyophilized powder). Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Native Chicken IgG
Synonyms:	Chicken Immunoglobulin G

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481