

Rabbit Anti-Avidin Secondary Antibody-Gold

Product Details

Ig Type:	IgG
Reactivity:	Avidin
Purification:	Protein A purified

Applications

Application:	ICA,IEM
Recommended	ICA: 1:20-200; IEM: 1:20-200

Properties

Stability & Storage:	Store at 2-8°C for 3-6 months, do not freeze.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Native Avidin
Synonyms:	Rabbit Anti-Avidin IgG H&L (Gold)

Research Background

The biological function of avidin is not known. It is a homotetrameric glycoprotein found in chicken egg white, and forms a strong non covalent specific complex with biotin (one molecule of biotin per subunit of avidin).

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481