

Rabbit Anti-Dog IgM Secondary Antibody

Product Details

Ig Type:	IgG
Reactivity:	Dog
Purification:	Protein A purified

Applications

Application:	ELISA,IF,IHC-F,IHC-P
Recommended	IHC-P: 1:100-500; IHC-F: 1:100-500; IF: 1:50-200; ELISA: 1:500-1000

Properties

Stability & Storage:	Store at -20°C or -80°C for 2 years (lyophilized powder). Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Native Dog IgM
Synonyms:	Dog Immunoglobulin M

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481