

## Rabbit Anti-Dog IgM Secondary Antibody-Gold

### Product Details

Ig Type: IgG  
Reactivity: Dog  
Purification: Protein A purified

### Applications

Application: ICA,IEM  
Recommended ICA: 1:20-200; IEM: 1:20-200

### Properties

Stability & Storage: Store at 2-8°C for 3-6 months, do not freeze.  
Shipping: Shipping with blue ice.

### Antigen Details

Immunogen: Native Dog IgM  
Synonyms: Immunoglobulin M;Rabbit Anti-Dog IgM (Gold)

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481