

## Rabbit Anti-Dog IgM Secondary Antibody-HRP

### Product Details

Ig Type: IgG  
Reactivity: Dog  
Purification: Protein A purified

### Applications

Application: ELISA,IHC-F,IHC-P,WB  
Recommended WB: 1:1000-10000; IHC-P: 1:100-500; IHC-F: 1:100-1000; ELISA: 1:1000-10000

### Properties

Stability & Storage: Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.  
Shipping: Shipping with blue ice.

### Antigen Details

Immunogen: Native Dog IgM  
Synonyms: Rabbit Anti-Dog IgM (HRP);Immunoglobulin M

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481