

Rabbit Anti-Dog IgM Secondary Antibody-PE-Cy7

Product Details

Ig Type:	IgG
Reactivity:	Dog
Purification:	Protein A purified

Applications

Application:	FCM,ICC/IF,IF
Recommended	IF: 1:100-1000; FCM: 1:100-1000; ICC/IF: 1:100-1000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Native Dog IgM
Synonyms:	Immunoglobulin M;Rabbit Anti-Dog IgM (PE-Cy7)

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481