

AffiniPure Goat Anti-Mouse IgG H&L (min X Hu, Bov, Hrs, Rb, Sw Sr Prot) Secondary Antibody-FITC

Product Details

Ig Type: IgG
Reactivity: Mouse
Purification: Protein G purified

Applications

Application: FCM,ICC/IF,IF
Recommended ICC/IF: 1:100-1000; IF: 1:100-1000; FCM: 1:100-1000

Properties

Stability & Storage: Store at -20°C or -80°C for 2 years (lyophilized powder). Avoid repeated freeze-thaw cycles.
Shipping: Shipping with blue ice.

Antigen Details

Immunogen: Native Mouse IgG
Synonyms: AffiniPure Goat Anti-Mouse IgG (H+L)(min X Hu, Bov, Hrs, Rb, Sw Sr Prot)

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481