

human IgG-Fc Secondary Antibody-Biotin (5I826)

Product Details

Ig Type:	IgG1
Reactivity:	Human
Clone:	5I826
Purification:	Protein A purified

Applications

Application:	ELISA,ICC/IF,IF,IHC-F,IHC-P,WB
Recommended	WB: 1:1000-20000; IHC-P: 1:100-1000; IHC-F: 1:100-1000; ELISA: 1:1000-20000; ICC/IF: 1:100-1000; IF: 1:100-1000; FCM: 1:100-1000

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	human IgG-Fc
Synonyms:	human IgG Fc;human IgG Fc fragment;Immunoglobulin G Fc fragment

Research Background

The Fc region is the constant region on an immunoglobulin molecule, the area that is exactly the same on all antibodies. The Fc region is found on the heavy chains and is not involved in binding antigens.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481