

Rat Anti-Mouse IgD Secondary Antibody (5R874)

Product Details

Ig Type:	IgG2a, κ
Reactivity:	Mouse
Clone:	5R874
Purification:	Protein A purified

Applications

Application:	FCM
Recommended	FCM: 1:100-500

Properties

Stability & Storage:	Store at 2°C-8°C, do not freeze.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Mouse IgD
Synonyms:	Mouse Immunoglobulin D

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481