


## 2'-Deoxy-N6-(4-methoxybenzyl)adenosine

## Chemical Properties

CAS No. :	157640-47-6
Formula:	C <sub>18</sub> H <sub>21</sub> N <sub>5</sub> O <sub>4</sub>
Molecular Weight:	371.39
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	2'-Deoxy-N6-(4-methoxybenzyl)adenosine is a Nucleoside Derivative - 6-Modified purine nucleoside.
Targets(IC50)	Nucleoside Antimetabolite/Analog,Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.6926 mL	13.4629 mL	26.9259 mL
5 mM	0.5385 mL	2.6926 mL	5.3852 mL
10 mM	0.2693 mL	1.3463 mL	2.6926 mL
50 mM	0.0539 mL	0.2693 mL	0.5385 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481