

L 735489

Chemical Properties


CAS No. : 147977-21-7

Formula: C₃₂H₅₁N₃O₆S

Molecular Weight: 605.83

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	L 735489 is a bioactive chemical.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.6506 mL	8.2531 mL	16.5063 mL
5 mM	0.3301 mL	1.6506 mL	3.3013 mL
10 mM	0.1651 mL	0.8253 mL	1.6506 mL
50 mM	0.033 mL	0.1651 mL	0.3301 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481