


## AMPA receptor antagonist-3

## Chemical Properties

CAS No. :	923272-18-8
Formula:	C <sub>20</sub> H <sub>19</sub> N <sub>5</sub> O <sub>2</sub> S
Molecular Weight:	393.46
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	AMPA receptor antagonist-3 is an AMPA receptor antagonist.
Targets(IC50)	iGluR

## Solubility Information

Solubility	DMSO: 30 mg/mL (76.25 mM), Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.5416 mL	12.7078 mL	25.4155 mL
5 mM	0.5083 mL	2.5416 mL	5.0831 mL
10 mM	0.2542 mL	1.2708 mL	2.5416 mL
50 mM	0.0508 mL	0.2542 mL	0.5083 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

## Reference

US20070027143A1

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481