


Sesquicannabigerol

Chemical Properties

CAS No. :	1334308-58-5
Formula:	C ₂₆ H ₄₀ O ₂
Molecular Weight:	384.60
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Sesquicannabigerol is classified as a phytocannabinoid.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.6001 mL	13.0005 mL	26.001 mL
5 mM	0.520 mL	2.6001 mL	5.2002 mL
10 mM	0.260 mL	1.3001 mL	2.6001 mL
50 mM	0.052 mL	0.260 mL	0.520 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481