

Phenmetrazine hydrochloride

Chemical Properties

CAS No. :	1707-14-8
Formula:	C11H16ClNO
Molecular Weight:	213.70
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	Phenmetrazine hydrochloride can induce self-administration behavior in dogs.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.6795 mL	23.3973 mL	46.7946 mL
5 mM	0.9359 mL	4.6795 mL	9.3589 mL
10 mM	0.4679 mL	2.3397 mL	4.6795 mL
50 mM	0.0936 mL	0.4679 mL	0.9359 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481