

YX0597

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Keep away from direct sunlight

Storage:

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	YX0597 is a potent and selective CDK9 PROTAC degrader that reduces RNA polymerase II S2 phosphorylation and MCL1 expression in GA0518 and AGS cells. It effectively inhibits the growth of GEAC cells (including radioresistant tumors) and suppresses tumor proliferation, metastasis, and cancer stem cells.
Targets(IC50)	CDK,PROTACs

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481