

LLC0424

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Keep away from direct sunlight
Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	LLC0424 is a potent and selective degrader of the nuclear receptor-binding SET domain 2 (NSD2), with a DC50 of 20 nM. It reduces H3K36me2 levels and inhibits the growth of acute lymphoblastic leukemia (ALL) cell lines, making it useful for research in ALL.
Targets(IC50)	Histone Methyltransferase,PROTACs

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481