


(3S)-3-Carboxy-3-hydroxypropanoyl-CoA

Chemical Properties

CAS No. : 52691-41-5
 Formula: C₂₅H₄₀N₇O₂₀P₃S
 Molecular Weight: 883.61
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	(3S)-3-Carboxy-3-hydroxypropanoyl-CoA [(3S)-3-Carboxy-3-hydroxypropanoyl-coenzyme A] is a type of 3-hydroxyacyl-CoA.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.1317 mL	5.6586 mL	11.3172 mL
5 mM	0.2263 mL	1.1317 mL	2.2634 mL
10 mM	0.1132 mL	0.5659 mL	1.1317 mL
50 mM	0.0226 mL	0.1132 mL	0.2263 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481