


Demycarosyl platenomycin

Chemical Properties

CAS No. :	52044-89-0
Formula:	C ₃₁ H ₅₁ N ₁ O ₁₁
Molecular Weight:	613.74
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Demycarosylplatenomycin (DM-PLM) is an antibiotic that exhibits mild activity against Gram-positive bacteria.
Targets(IC50)	Antibacterial, Antibiotic

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.6294 mL	8.1468 mL	16.2935 mL
5 mM	0.3259 mL	1.6294 mL	3.2587 mL
10 mM	0.1629 mL	0.8147 mL	1.6294 mL
50 mM	0.0326 mL	0.1629 mL	0.3259 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481