

## Cy3 Dextran (MW 40000)

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: **Keep away from direct sunlight**  
**Powder: -20°C for 3 years | In solvent: -80°C for 1 year**  
Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Cy3Dextran (MW 40000) is a fluorescent labeling agent that combines Cy3 fluorescent dye with Dextran. The Cy3 fluorophore is commonly utilized in applications such as immunolabeling, nucleic acid labeling, fluorescence microscopy, and flow cytometry. Dextran serves as a plasma volume expander and inhibits platelet aggregation and coagulation factors [Ex/Em = 550/570 nm].
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481