

m-PEG-PEI800-Cy5

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Keep away from direct sunlight
Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	m-PEG-PEI800-Cy5 is a fluorescent dye comprising CY5, methoxy polyethylene glycol (mPEG), and polyethyleneimine (PEI). This compound is utilized in gene delivery, nano-drug delivery, cellular imaging, biosensing, and bioanalysis applications, with excitation/emission (Ex/Em) wavelengths of 633/670 nm.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481