

HO-PEG-SCM (MW 3400)

Chemical Properties

CAS No. : 196936-11-5

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	HO-PEG-SCM (MW 3400) is a heterobifunctional hydroxy PEG succinimidyl carboxymethyl ester, commonly used as a crosslinking agent or spacer between two distinct chemical entities.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481