

HPK1-IN-58

## Chemical Properties

CAS No. :

Formula: C<sub>26</sub>H<sub>30</sub>N<sub>8</sub>O

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	HPK1-IN-58 is an inhibitor of HPK1 (IC <sub>50</sub> : 2.6 nM) and SLP76 (IC <sub>50</sub> : 20 nM). It enhances IL-2 secretion, thereby reversing PGE <sub>2</sub> -induced immunosuppression. HPK1-IN-58 is applicable for research in antitumor immunity.
Targets(IC <sub>50</sub> )	MAPK,Interleukin

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481