

## (3S)-Hydroxytetradecanedioyl-CoA

### Chemical Properties

CAS No. :

Formula: C<sub>35</sub>H<sub>60</sub>N<sub>7</sub>O<sub>20</sub>P<sub>3</sub>S

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	(3S)-Hydroxytetradecanedioyl-CoA, also known as (3S)-Hydroxytetradecanedioyl-coenzyme A, is a specific type of (S)-3-hydroxyacyl-coenzyme A.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481