


(S)-3-Hydroxyhexacosanoyl-CoA

Chemical Properties

CAS No. :	1245945-31-6
Formula:	C47H86N7O18P3S
Molecular Weight:	1162.21
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	(S)-3-Hydroxyhexacosanoyl-CoA ((S)-3-Hydroxyhexacosanoyl-coenzyme A) is a type of very long-chain hydroxyacyl-CoA.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.8604 mL	4.3021 mL	8.6043 mL
5 mM	0.1721 mL	0.8604 mL	1.7209 mL
10 mM	0.086 mL	0.4302 mL	0.8604 mL
50 mM	0.0172 mL	0.086 mL	0.1721 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481