

2-Benzoylsuccinyl-CoA

Chemical Properties


CAS No. : 1383686-27-8

Formula: C32H44N7O20P3S

Molecular Weight: 971.71

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	2-Benzoylsuccinyl-CoA (2-Benzoylsuccinyl-coenzyme A) is a type of acyl-coenzyme A.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.0291 mL	5.1456 mL	10.2911 mL
5 mM	0.2058 mL	1.0291 mL	2.0582 mL
10 mM	0.1029 mL	0.5146 mL	1.0291 mL
50 mM	0.0206 mL	0.1029 mL	0.2058 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481