

## DOT1L ligand-1

## Chemical Properties

CAS No. :

Formula: C<sub>24</sub>H<sub>19</sub>ClF<sub>2</sub>N<sub>8</sub>O<sub>7</sub>S

Molecular Weight:

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	DOT1Lligand-1 is a ligand for DOT1L. It serves as a ligand for target proteins (Ligands for Target Protein for PROTAC) and is utilized in the synthesis of PROTACs targeting DOT1L, such as MS2133.
Targets(IC50)	Others,Ligands for Target Protein for PROTAC

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481