

Cy7.5 HA MA (MW 100000)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Keep away from direct sunlight
Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Cy7.5HA MA (MW 100000) is a methacrylate derivative of hyaluronic acid labeled with CY7.5, extensively utilized in the fields of bio-labeling and imaging, drug delivery, as well as tissue engineering and regenerative medicine (Ex/Em = 770/820 nm).
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481