

Anti-STEAP2 Antibody

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: -20°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Anti-STEAP2 Antibody (AZD0516 antibody) is a selective inhibitor of the prostate-specific six transmembrane epithelial antigen 2 (STEAP2). It inhibits the proliferation, migration, and invasion of prostate cancer cells by blocking the STEAP2-mediated ERK/MAPK signaling pathway. Anti-STEAP2 Antibody shows potential for prostate cancer research.
Targets(IC50)	ERK,p38 MAPK

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481