


## 2-Propylpent-4-enoyl-CoA

## Chemical Properties

CAS No. :	148963-27-3
Formula:	C <sub>29</sub> H <sub>48</sub> N <sub>7</sub> O <sub>17</sub> P <sub>3</sub> S
Molecular Weight:	891.71
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	2-Propylpent-4-enoyl-CoA (2-Propylpent-4-enoyl-coenzyme A) is a derivative of coenzyme A.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.1214 mL	5.6072 mL	11.2144 mL
5 mM	0.2243 mL	1.1214 mL	2.2429 mL
10 mM	0.1121 mL	0.5607 mL	1.1214 mL
50 mM	0.0224 mL	0.1121 mL	0.2243 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481