


(S)-Carnitiny-CoA, betaine

Chemical Properties

CAS No. : 1202823-44-6
 Formula: C₂₈H₅₀N₈O₁₈P₃S
 Molecular Weight: 911.73
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	(S)-Carnitiny-CoA, betaine ((S)-Carnitiny-coenzyme A, betaine) is a type of 3-hydroxyacyl-coenzyme A.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.0968 mL	5.4841 mL	10.9682 mL
5 mM	0.2194 mL	1.0968 mL	2.1936 mL
10 mM	0.1097 mL	0.5484 mL	1.0968 mL
50 mM	0.0219 mL	0.1097 mL	0.2194 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481