

(S)-E3 Ligase Ligand-linker Conjugate 140

Chemical Properties

CAS No. :

Formula: C₃₀H₃₄N₆O₅

Molecular Weight:

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	(S)-E3LigaseLigand-linker Conjugate 140 is an E3 ubiquitin ligase ligand-linker conjugate (E3ligaseligand-linker conjugate) that is applicable in the synthesis of GDC-2992.
Targets(IC50)	E3 Ligase Ligand-Linker Conjugates

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481